




<b>Product</b> <b>Item number</b>	VALUE folding shower seat with supporting leg, fixed height R1602000
<b>Color</b>	VVS:777526435 NRF:6191657 RSK:7973958 EAN:5708590185927 HMI: 25572 NBS:N13:448 White
<b>Product description</b>	 <p>VALUE folding shower seat and backrest create a simple solution for users who need to sit down while showering.</p> <p>Folding shower seat with support leg, seat pale gray Incl. wall mounting kit</p>
<b>Adjustability</b>	Fixed mounting
<b>Product capacity</b>	Max. vertical load: 275 lbs.
<b>Materials</b>	Diecast, powder-coated aluminium. Rilsan-coated steel, which is comfortable to the touch, easy to clean, and protects the metal against corrosion. Plastic part is made of vacuum-moulded polyethylene (PE). <i>Flammability-tested (ISO 8191-1 and ISO 8191-2). Test report available.</i> Toilet seat and cover in thermoset plastic (urea).
<b>Ergonomy</b>	Ergonomic seat with a comfortable surface.
<b>Features</b>	The seat can fold up when not in use or to create more space. Leg for extra support. With run-off grooves.
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



# PRESSALIT®

<b>Product</b>	VALUE folding shower seat with supporting leg, fixed height
<b>Item number</b>	R1602000
<b>Cleaning</b>	<p>This product should be cleaned with hot water or a mild cleaning solution. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p><b>DISINFECTANTS:</b> CaviCide</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>
<b>Declaration of conformity</b>	<p>This product is in compliance with all relevant requirements of regulation (EU) 2017/745 on medical devices.</p>

[pressalit.com](http://pressalit.com)

## PRESSALIT®