



|   |  |
|---|--|
| <b>Product</b><br><b>Item number</b>            | CareSystem shower seat 310, manually height adjustable<br>R9705000   |
| <b>Colour</b>                                   | White  |
| <b>Product description</b>                      | Shower seat 310, "soft feel", height adjustable, without vertical track.<br>Discontinued as per 31/12-2015<br><br>Shower seat 310, "soft feel", height adjustable, without vertical track.<br>Discontinued as per 31/12-2015   |
| <b>Adjustability</b><br><b>Product capacity</b> | Height adjustable<br><br>Max. vertical load: 150 kg.   |
| <b>Materials</b>                                | Basic material: 20 micrometre anodized aluminium<br><br>Seat: Polyurethane foam with marine ply core<br><br>Seat surface: Gum lac<br><br>Lever: Stainless steel w/nylon knob   |
| <b>Ergonomy</b>                                 | Soft and slip-resistant material with a warm feel.<br><br>Ergonomically designed seat.   |
| <b>Features</b>                                 | The seat can fold up when not in use or to make more space.<br><br>Height adjustable.<br><br>Without vertical track.   |
| <b>Cleaning</b>                                 | Pressalit recommends that chalky deposits be removed from wall tracks and other aluminum parts using a solution of three parts of 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the parts should be wiped over with a damp cloth wrung out in clear water. |



# PRESSALIT®

|                    |  |
|--------------------|--|
| <b>Product</b>     | CareSystem shower seat 310, manually height adjustable |
| <b>Item number</b> | R9705000   |

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®