




PRESSALIT®

Product	Toilet paper holder for TMA fold down grab bars		
Item number	R1175000		
Color	VVS:615723950 NRF:6191655 EAN:5708590211725 HMI: 35271 NBS:N13:436		
	White		
Product description	<p>Paper holder for TMA arm supports. Fixed, can be mounted on the left or right handle.</p> <p>Paper holder for TMA arm supports. Fixed, can be mounted on the left or right handle.</p>		
Adjustability	Fixed mounting		
Materials	Powder-coated stainless steel.		
Ergonomy	Can be mounted at the right or the left handle.		
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.		

pressalit.com

PRESSALIT®