




# PRESSALIT®

<b>Product</b>	Wall-mounted toilet 21.3"
<b>Item number</b>	R2091000
<b>Color</b>	VVS:613230000 NRF:6190522 EAN:5708590213897 HMI: 35273 NBS:N13:10,11,300 White
<b>Product description</b>	 <p>Toilet, wall-mounted, distance from wall to the front: 21.3" Excl. mounting kit. For mounting directly on wall, a mounting kit from V&amp;B must be ordered from Pressalit at an additional cost. For mounting on built-in cistern or Pressalit toilet lifter, mounting kit is included in the delivery.</p> <p>Toilet, wall-mounted, distance from wall to the front: 21.3". Excl. mounting kit</p>
<b>Product capacity</b>	Max load: 880 lbs.
<b>Materials</b>	Ceramic
<b>Ergonomy</b>	Wall-mounted toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a floor lift, and the caregiver's feet.
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.

[pressalit.com](http://pressalit.com)

## PRESSALIT®