




PRESSALIT®

Product Item number	Modular shelves for PLUS sink bracket, 23.6" rod RT481
Product description	<p>VVS:777558260 NRF:7501026 RSK:7973954 EAN:5708590146928 HMI: 124127 NBS:N13:449</p>  <p>PLUS modular shelves are practical shelves attached to a pole, they are available for fitting to PLUS wall tracks or PLUS wash basin brackets, as well as for wall-mounting.</p> <p>Modular shelf system: Rod 23.6", for mounting on bracket for sink, incl. two shelves 15.8", towel rail and two trays.</p>
Adjustability Product capacity	<p>Height adjustable</p> <p>Maximum load: 24 lbs. Load per shelf: 8 lbs. Load on towel holder: 8 lbs.</p>
Materials Features Warranty	<p>Aluminium. Stainless steel. Laminate.</p> <p>Height adjustable.</p> <p>Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.</p>



PRESSALIT®

Product	Modular shelves for PLUS sink bracket, 23.6" rod
Item number	RT481

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®