



# PRESSALIT®

<b>Product</b>	Floor standing toilet with S-trap, extra high
<b>Item number</b>	R2097000
	VVS:605034414 NRF:6190523 EAN:5708590335346 HMI: 105699
<b>Color</b>	White
<b>Product description</b>	<p>Toilet, floor-standing, extra high 18.1". With concealed S-trap and dual flushing system preset at 3/6 liters. Can be changed to 2/4 liters Excl. mounting kit</p> <p>Toilet, floor-standing, extra high 18.1". With concealed S-trap and dual flushing system preset at 3/6 liters. Can be changed to 2/4 liters Excl. mounting kit</p>
<b>Materials</b>	Ceramic
<b>Ergonomy</b>	A high toilet can with advantage be used by persons who would otherwise need a raised toilet seat.
<b>Features</b>	Extra high. With concealed S-trap and dual flushing system
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.

[pressalit.com](http://pressalit.com)

# PRESSALIT®